



*Fig. 1 Proximal humerus locking plate failure with joint penetration.*



*Fig. 2 Locked construct failure with new tibia shaft fracture.*



*Fig. 3 Bridge plating with plate fatigue failure through unplugged locking screw hole.*